



Product datasheet



Hypred quick spray

Sprayable liquid for teat disinfection.





Advantages of HYPRED QUICK SPRAY

- Very powerful disinfectant with LSA®
- Disinfection of the teat and udder base
 - Gentle cleaning of the teats



Teat protection

Hypred Quick Spray is enriched on hydrating and emollient substances bringing exceptional cosmetic virtues.

Teat disinfection

Hypred Quick Spray is a disinfectant based on LSA® (hydroxy acids complex).

The LSA® complex is effective in 30 seconds on the germs that cause mastitis: Staphylococcus aureus, Streptococcus uberis, Streptococcus Agalactiae, Escherichia coli, ... (EN 1656).

The LSA® complex is also effective on Orthopoxvirus responsible of the vaccinia (EN 14675) and on the yeasts Candida albicans (EN 1657).



Direction for use

- Ready to use before and after milking
- -Before milking pulverize Hypred Quick Spray on the teat, let it act, and wipe with a paper of single use or with a clean cloth.
- Aftermilking apply Hypred Quick Spray.

Main components of Hypred quick spray

-LSA® (complex of hydroxy-acid), hydrating, softening, healing, tensoactives

- Active(s) substance(s) per 100g of product: acid L - (+) - lactic 2g + salicylic acid 0.099g

Disinfectants are regulatory biocides. They have guarantees of efficiency and protection of human, animals and the environment.

Use biocides safely. Before use read the label and product information.



